

Reverse Transfection of Arrayed crRNA

This protocol is for reverse transfection of Dharmacon Edit-R crRNA arrayed libraries.

Table of Contents

1. Reverse transfection of arrayed crRNAs	1
Product description	1
Materials.....	1
Guidelines for resuspension of the crRNA	2
Reverse transfection of arrayed crRNA.....	2
2. Appendix	3
Optimization of transfection conditions for delivery of crRNA:tracrRNA.....	3
Gene editing assay recommendations.....	3
3. Frequently asked questions	4

1 Reverse transfection of arrayed crRNAs

Product description

Libraries of predesigned synthetic crRNAs provide an opportunity to apply the CRISPR-Cas9 system for functional gene knockout analysis in an arrayed format. Dharmacon™ Edit-R™ crRNA libraries consist of synthetic crRNAs that are predesigned using the proprietary Edit-R CRISPR RNA algorithm. This algorithm was trained on functional knockout data and has demonstrated the ability to select guide RNA target regions more likely to give functional knockout of the protein, not just create a DNA cut. The algorithm additionally includes specificity scoring using an internal alignment tool for complete off-target identification, identifying both mismatches and gaps.

Edit-R crRNA Libraries are available as the following:

- Catalog libraries of predefined gene family collections for human and mouse
 - » Four crRNA per gene, 80 wells per plate, columns 1 and 12 left empty, 0.5 nmol per well in 96-well plates
 - » Provided in NUNC Polystyrene 96-well V-bottom plates (Cat #249952)
- Cherry-pick crRNA libraries are created by a researcher from their own gene list. Learn more or get started at <http://dharmacon.gelifsciences.com/cherry-pick-libraries/>
 - » One to five crRNA per gene (minimum of three is recommended) in 96-well plates, minimum 40 wells
 - » Customizable plate layout; Edit-R catalog control crRNAs may be added to any wells within the plate(s).
 - » Provided in NUNC Polystyrene 96-well V-bottom plates (Cat #249952)

Materials

- Plates of synthetic crRNAs, up to 2 nmol per well in 96-well plates

Additional required materials are listed below and are not provided with crRNA library purchase:

Edit-R tracrRNA, 5, 20 or 50 nmol (Dharmacon, Cat #U-002005-XX). tracrRNA is required for use with all synthetic crRNA reagents.

- DharmaFECT transfection reagent (formulation is dependent on specific cell line of interest)
 - » DharmaFECT 1 Cat #T-2001-XX
 - » DharmaFECT 2 Cat #T-2002-XX
 - » DharmaFECT 3 Cat #T-2003-XX
 - » DharmaFECT 4 Cat #T-2004-XX
- 10 mM Tris pH 7.4 nuclease-free buffer solution (Dharmacon, Cat #B-006000-100)
- 96-well tissue culture plates
- 96-well V-bottom polystyrene standard storage plates or deep well plates (for example, NUNC Cat #249952 or Cat #12-565-553)
- Assay for cell viability such as CellTiter-Blue® Cell Viability Assay (Promega Corp., Cat #G8081)
- Assay(s) for detecting gene editing events in a cell population (<http://dharmacon.gelifsciences.com/uploadedFiles/Resources/crrna-positive-controls-protocol.pdf>)
- Positive control crRNA and detection primers for assessment of gene editing (<http://dharmacon.gelifsciences.com/gene-editing/crispr-cas9/edit-r-synthetic-positive-crrna-controls-and-detection-primers/>)
- Non-targeting control crRNA (<http://dharmacon.gelifsciences.com/gene-editing/crispr-cas9/edit-r-synthetic-crrna-non-targeting-controls/>)
- Assay-specific positive control crRNA (defined by researcher)
- Growth medium: antibiotic-free cell culture medium (with serum and/or supplements) recommended for maintenance of the cells of interest
- Serum-free and antibiotic-free cell culture medium for preparation of transfection mix (for example, MEM-RS, HyClone Cat # SH30564).



For phenotypic analysis with the arrayed synthetic crRNA libraries we strongly recommend using cell lines that constitutively express Cas9 nuclease. Transfection of crRNA:tracrRNA into a cell line that is constitutively expressing Cas9 nuclease results in a higher percentage of edited cells thus allowing for easier downstream high-throughput phenotypic analysis. For generation of the Cas9 stable expressing cells please follow recommendations in the [Gene Engineering with Lentiviral Cas9 Particles and Synthetic CRISPR RNAs](#) manual.



The phenotypic analysis, including assay optimization and analysis, is cell line and assay specific and requires optimization by the researcher.

Guidelines for resuspension of the crRNA

1. The crRNA libraries are shipped at ambient temperature as dry pellets of RNA in each well and should be stored at -20 °C upon arrival in a manual defrost or non-cycling freezer. If necessary, crRNAs can be stored as dry pellets (unopened) at 4 °C for several weeks.
2. Briefly centrifuge plates to ensure that the crRNA is collected at the bottom of the well.
3. Wipe adhesive foil cover with 70% ethanol or other RNase-decontamination solution.
4. Pierce or carefully peel back the foil seal to gain access to wells. Use caution and avoid shredding the seal.
5. If you are starting with a plate of 0.5 nmol per well, resuspend arrayed crRNAs to 10 µM solution by adding 50 µL nuclease-free 10 mM Tris pH 7.4 buffer to 0.5 nmol of crRNA (for different quantities of crRNAs see Table 3).
6. Pipette solution up and down 3-5 times while avoiding introduction of bubbles.
7. Seal the plate and place it on an orbital mixer/shaker for 70-90 minutes at room temperature.
8. Briefly centrifuge plates to collect solution to bottom of the wells.
9. From this master crRNA plate, generate plates with 1 µM working concentration of crRNAs using nuclease-free, 10 mM Tris pH 7.4 buffer. This eliminates the subsequent requirement for pipetting of very small volumes.
10. crRNA plates may now be used immediately, aliquoted into single-use plates or stored at -20 °C in a manual defrost or non-cycling freezer.
11. For storage, seal plates with appropriate adhesive or heat seals.

Reverse transfection of arrayed crRNA

The following is a general protocol for reverse transfection of arrayed crRNA libraries using stable Cas9-expressing mammalian cells. For guidelines using arrayed crRNA plates with forward transfection see <http://dharmacon.gelifsciences.com/uploadedFiles/Resources/edit-r-crrna-libraries-guide-protocol.pdf>. Optimal plating density will depend on growth characteristics of specific cell lines and assay requirements and these parameters should be determined experimentally. Exact parameters for crRNA:tracrRNA transfection in your cells of interest should be empirically determined through careful optimization prior to experimentation [see Appendix for Optimization of transfection conditions with crRNA:tracrRNA](#). Catalog crRNA library plates are supplied with columns 1 and 12 empty to allow addition of researcher-defined controls. We suggest including the following controls:

- Untreated cells
- Positive control crRNA
- Negative control: non-targeting crRNA

Day 1

This protocol is for transfection of one arrayed crRNA 96-well plate, in triplicate, at final 25 nM concentration of the crRNA:tracrRNA. Calculations are done for quadruplicates providing excess for the ease of pipetting.

1. Transfer 10 µL of 1 µM working crRNA solution to each well of a 96-well V-bottom transfection mixture plate (NUNC polystyrene 96 well V- bottom plates Cat #249952 or other appropriate plates may be used).



Positive and negative crRNA controls can be added to empty wells of the V-bottom transfection mixture plate (columns 1 and 12 in catalog libraries).

2. Resuspend Edit-R tracrRNA to 10 µM stock solution in nuclease-free, 10 mM Tris pH 7.4 buffer (for example, add 500 µL to 5 nmol of tracrRNA).
3. Prepare 333 nM tracrRNA working solution by adding 120 µL of 10 µM tracrRNA stock solution to serum-free medium for a total volume of 3.6 mL. This will allow preparation of one arrayed crRNA library plate in triplicate (total of three plates) and includes excess for ease of pipetting.
4. Add 30 µL of 333 nM working tracrRNA solution to each well of 96-well V-bottom transfection mixture plate containing 10 µL of crRNA, prepared in step 1. This will bring the concentration of the crRNA and tracrRNA to 250 nM.
5. Prepare transfection reagent working solution by diluting the transfection reagent stock solution in serum-free medium for a total volume of 5 mL. This volume will allow preparation of one crRNA library plate in triplicate and includes excess for ease of pipetting. For example, if the optimal amount of transfection reagent was determined to be 0.05 µL per well of cells, add 25 µL of transfection reagent stock solution to serum-free medium for a total volume of 5 mL. For preparations of other transfection reagent concentrations see Table 1.
6. Add 40 µL of transfection reagent working solution to each well of 96-well V-bottom transfection mixture plate containing the crRNA:tracrRNA complex. This brings the total volume to 80 µL.
7. Immediately mix by pipetting gently up and down and incubate for 20 minutes at room temperature.
8. During this incubation time, cells may be prepared for seeding. Trypsinize cells, spin down, remove medium and resuspend in cell growth medium such that when 80 µL are added per well it will result in the previously determined optimal cell density for your transfection. For example, if optimal cell density was determined to be 10,000 cells per well, prepare 28 mL of cell suspension at 1.25×10^5 cells/mL for three plates.
9. Briefly mix (by pipetting) the transfection mix prepared in step 7 and add 20 µL of the transfection mixture to each well of three 96-well tissue culture plates.
10. Add 80 µL of the cell suspension prepared in step 8 to each well of the three 96-well tissue culture plates. This will bring the volume to 100 µL and the final concentration of the crRNA:tracrRNA complex to 25nM.
11. Incubate transfected cells at 37 °C in a humidified CO₂ incubator for 48-72 hours before proceeding with the phenotypic assay or gene editing analysis.

Table 1. Transfection reagent amounts for transfection of one 96-well crRNA library plate in triplicate.

Transfection reagent volume per well of cells (μL)	Transfection reagent volume (μL)	Serum-free medium volume (mL)
0.01	5	5.00
0.025	12.5	4.99
0.05	25	4.98
0.1	50	4.95
0.2	100	4.90

2 Appendix

Optimization of transfection conditions for delivery of crRNA:tracrRNA

To obtain the highest transfection efficiency of the Edit-R crRNA:tracrRNA components with minimal effects on cell viability, we recommend carefully optimizing transfection conditions for each cell line using a positive control crRNA. The transfection optimization can be easily performed in a 96-well format allowing for testing of multiple transfection conditions. Transfection conditions that have previously been optimized for siRNA delivery are a reasonable starting point for crRNA:tracrRNA transfection optimization. Cell type-specific guidelines for the four DharmaFECT transfection reagent formulations can be found in the [DharmaFECT Cell Type Guide](#).

 **The transfection reagent amounts for reverse transfection are typically different from forward transfection, and optimal conditions should be empirically determined for your cells.**

The optimization experiment should include two to three cell densities and a range of DharmaFECT transfection reagent volumes. Our recommendations for the different components in the transfection optimization experiment are as follows:

- 0.01 to 0.4 μL/well of DharmaFECT 1, 2, 3, or 4 transfection reagent in a 96-well plate
- 25 nM positive control (PPIB or DNMT3B) crRNA:tracrRNA per well (recommended range 20 nM-50 nM).

Use Table 2 for guidance in preparation of samples for 96-well transfection optimization.

Table 2. Preparing samples for 96-well transfection optimization. Mix the diluted crRNA:tracrRNA and diluted transfection reagent and incubate the transfection mixture for 20 minutes. Transfer the transfection mixture to a well of a tissue culture plate and add 80 μL of cell suspension in growth medium at the desired cell density. This will bring the total volume per well to 100 μL and the concentration of crRNA:tracrRNA complex to 25 nM. Volumes (μL) are shown per ONE well of a 96-well plate; for triplicate wells multiply all values by 3.5 to have sufficient volume for three wells and to account for pipetting error. For the diluted transfection reagent prepare a larger volume to enable accurate pipetting of the small volumes that are required per one well.

Transfection condition	Tube 1: Diluted crRNA:tracrRNA			Tube 2: Diluted DharmaFECT transfection reagent		Final reagent volumes		
	Serum-free medium (μL)	Volume of 1 μM crRNA (μL)	Volume of 1 μM tracrRNA (μL)	Serum-free medium (μL)	Volume of DharmaFECT transfection reagent (μL)	Transfection mixture volume (μL)	Cell suspension in Growth medium (μL)	Total volume per 96-well (μL)
0.01 μL/well	5	2.5	2.5	9.99	0.01	20	80	100
0.025 μL/well	5	2.5	2.5	9.975	0.025	20	80	100
0.05 μL/well	5	2.5	2.5	9.95	0.05	20	80	100
0.1 μL/well	5	2.5	2.5	9.9	0.1	20	80	100
0.2 μL/well	5	2.5	2.5	9.8	0.2	20	80	100
0.3 μL/well	5	2.5	2.5	9.7	0.3	20	80	100
0.4 μL/well	5	2.5	2.5	9.6	0.4	20	80	100
Untreated (0 μL/well)	10	0	0	10	0	20	80	100

At 48-72 hours post-transfection, perform a cell viability assay to determine the highest lipid concentration that has minimal cell toxicity (≥ 80% cell viability is preferred). After assaying for cell viability, we recommend that you carefully wash the cells once with PBS and proceed with a gene editing assay (see below) to determine the condition that produces the best editing efficiency. Use these optimal conditions for subsequent transfection of your selected Cas9-expressing cell lines with the Edit-R crRNA:tracrRNA.

Gene editing assay recommendations

A commonly used method for detection of insertions and deletions (indels) in a cell population is a DNA mismatch detection assay such as using T7 Endonuclease I (T7EI). The mismatch detection assay can be performed on either purified genomic DNA or whole cell lysate. For a detailed protocol see <http://dharmacon.gelifsciences.com/uploadedFiles/Resources/crrna-positive-controls-protocol.pdf>.

3 Frequently asked questions

How should I store my crRNA and tracrRNAs?

RNA oligonucleotides should be stored at -20 °C or -80 °C in a non-frost free freezer, either as a dried pellet or resuspended in an RNase-free 10 mM Tris solution buffered to pH 7.4 to help with stability during freeze-thaw cycles. We recommend that RNA oligonucleotides be resuspended to a convenient stock concentration (Table 3) and stored in small aliquots to avoid multiple freeze-thaw cycles. RNA oligonucleotides should not go through more than four to five freeze thaw cycles. If degradation is a concern, the integrity of the RNA oligonucleotides can be evaluated on an analytical PAGE gel.

Table 3. Making stock solutions of crRNA and tracrRNA.

crRNA or tracrRNA amount (nmol)	Volume (µL) of 10 mM Tris pH 7.4 to be added for desired final concentration	
	100 µM stock	10 µM stock
0.5	-	50
2	20	200
5	50	500
20	200	2000
50	500	Exceeds tube volume*

*When tube volume is exceeded, make a 100 µM stock and dilute it 1:10 to obtain the 10 µM stock.

What is the stability of the Edit-R crRNA and tracrRNAs?

Dried RNA oligonucleotide pellets are stable at room temperature for two to four weeks, but should be placed at -20 °C or -80 °C for long-term storage. Under these conditions, the dried crRNA and tracrRNAs will be stable for at least one year. Once resuspended the RNAs are stable for at least one year when stored at -20 °C or -80 °C in a non-frost free freezer. Maintaining sterile, RNase- and DNase-free conditions is always recommended as a critical precaution.

Can I use my siRNA transfection protocols to transfect Edit-R synthetic crRNA and tracrRNAs?

Previously optimized protocols to transfect synthetic siRNA into your cells of interest can be a good starting point for transfection of synthetic crRNA:tracrRNA.

Can I use a different transfection reagent other than DharmaFECT transfection reagents to deliver the Edit-R components into my cells?

We cannot predict the performance of other transfection reagents, nor can we troubleshoot experiments performed with any reagent other than DharmaFECT transfection reagents. However, other suitable transfection reagents designed for RNA transfection could be utilized provided transfection conditions are carefully optimized for each cell line of interest.

Can I co-transfect arrayed synthetic crRNA:tracrRNA with the Edit-R Cas9 Nuclease Expression plasmids?

You can perform genome engineering by transient transfection of the synthetic crRNA:tracrRNA with the Edit-R Cas9 Nuclease Expression plasmids using DharmaFECT Duo transfection reagent. However, for performing phenotypic analysis in the cell population in a high-throughput manner, we have found that crRNA:tracrRNA transfection in a cell line that stably expresses Cas9 nuclease produces higher efficiency gene editing in the cell population with lower toxicity associated with the transfection.

CellTiter-Blue® Cell Viability Assay is a trademark of Promega Corp. GE, imagination at work and GE monogram are trademarks of General Electric Company. Dharmacon is a trademark of GE Healthcare companies. All other trademarks are the property of General Electric Company or one of its subsidiaries. ©2017 General Electric Company—All rights reserved. Version published January 2016. GE Healthcare UK Limited, Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK



Orders can be placed at:
gelifesciences.com/dharmacon

Customer Support: cs.dharmacon@ge.com
Technical Support: ts.dharmacon@ge.com or
+1 800 235 9880; +1 303 604 9499 if you have any questions.