# Dharmacon™

RNAi, Gene Expression & Gene Editing

# Yeast Tet-promoters Hughes (yTHC) Parental Strain – R1158

Cat. #YSC1210

The Yeast Tet-promoters Hughes (yTHC) Parental Strain contains 632 essential yeast genes for which expression is regulated by doxycycline. The endogenous promoter has been replaced with a TET-titratable promoter in the genome. Thus, the expression of the gene can be switched off by the addition of doxycycline to the yeast's growth medium.

The yTHC mutant strains are provided in the haploid MATa strain R1158. This strain was created from the background strain BY4741 by a one-step integration of the tTA transactivator, under the control of the CMV promoter, at the URA3 locus (*Hughes et al., Cell.* 2000 Jul 7;102(1):109-26).

#### Storage:

- 4 °C for up to one week
- -80 °C indefinitely

#### **Product Description:**

Culture of S. cerevisiae in YPD broth + 15% glycerol.

### Making a Stock Culture:

Once the strain has been streaked, we recommend making a stock of the pure culture. Inoculate the pure culture in YPD broth and incubate for 48 hours at 30 °C. Transfer 850  $\mu$ L of culture into a sterile polypropylene tube and add 150  $\mu$ L sterile glycerol to make a 15% glycerol freezing solution. Vortex the culture to evenly mix the glycerol throughout the culture. The culture can be stored indefinitely at -80 °C.

## Genotype of S. cerevisiae R1158:

URA::CMV-tTA MATa his3-1 leu2-0met15-0l

GE, imagination at work and GE monogram are trademarks of General Electric Company. Dharmacon is a trademark GE Healthcare companies. All other trademarks are the property of General Electric Company or one of its subsidiaries. ©2014 General Electric Company—All rights reserved. Version published June 2014. GE Healthcare UK Limited, Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

