

siRNA Screening: Development of Hit Stratification Strategies


Žaklina Strezoska, Annaleen Vermeulen, Emily M. Anderson, Anja Smith, Devin Leake
Dharmacon, now part of GE Healthcare, 2650 Crescent Drive, Suite #100, Lafayette, CO 80026, USA

Introduction

While synthetic siRNA libraries are powerful tools for functional genomic screens, off-target effects mediated by siRNA seed interactions with the 3' UTR of unintended targets can result in false positives. Given the frequency of off-target effects in some assays, the development of hit validation/stratification strategies is imperative. In the following study we have compared two strategies for identification of high confidence hits: 1) a multiple reagent approach where two or more individual siRNAs induce the same phenotype and 2) a chemical modification approach where hit confirmation is achieved using pools of siRNA that contain specificity enhancing modifications. A comparison of these two strategies (using a collection of primary hits generated from a cell viability screen) reveals significant overlap between the high confidence hits identified. However, for low confidence hits, i.e. where a single siRNA induces a phenotype, the concern is that an important hit will be missed. To determine if the phenotype is due to gene targeting or a seed-mediated off-target effect, a chimeric approach was used whereby a gene-specific seed sequence is introduced into a non-targeting siRNA scaffold. Together, these data provide well-defined approaches for prioritization of hits derived from RNAi screens.


siRNA screening workflow

RNAi screens typically result in a large number of hits. Distinguishing between true hits and false positives due to off targets is a crucial task in first line hit validation.



Mechanism of RNAi-mediated effects

RNAi can lead to both specific target mRNA down-regulation and nonspecific off-target effects due to partial complementarity with unintended mRNAs through seed region matches to the 3' UTR. Typical results of microarray expression analysis is shown below.




Strategies for reducing off-targeting

Bioinformatics – Select for potent siRNAs and assess for seed content to promote specificity (1,4,5)


Chemical modifications – Expand sequence space enabling potent siRNA sequence selection and interfere with off-target mechanism (2)

Pooling – Promote competition between potent siRNAs for optimal silencing and dilution of off-targets associated with individual siRNAs (3)

Specificity-enhanced pools are ON-TARGETplus™, SMARTpool™, siRNA reagents, standard siRNA and pools are siGENOME SMARTpool and individual siRNA reagents




Goal of this study



Standard pool = siGENOME SMARTpool reagent (pool of 4 unmodified siRNAs)
Standard pool deconvolution = 4 individual siGENOME siRNAs that make up the pool
Specificity-enhanced pool = ON-TARGETplus SMARTpool reagent (pool of 4 modified siRNAs)

Standard and specificity-enhanced pools have comparable silencing



The target mRNA knockdown was compared between the standard and specificity-enhanced siRNA pools. Statistically significant silencing differences, $p < 0.01$ (*Specificity enhanced pool better, # Standard pool better)

Deconvolution of the standard siRNA pool vs. specificity enhanced pool

Standard pool	≥ 2 Standard siRNAs	Specificity enhanced pool	Confirmation
CENPE	CENPE	CENPE	Confirmed by both methods
COPB2	COPB2	COPB2	
MAPK3	MAPK3	MAPK3	Confirmed by one method
PDGFRB	PDGFRB	PDGFRB	
PIK4CA	PIK4CA	PIK4CA	Confirmed by neither method
STK6	STK6	STK6	
TESK1	TESK1	TESK1	Low or not expressed genes
FN3KRP	FN3KRP	FN3KRP	
CDC20		CDC20	Confirmed by one method
CDK11		CDK11	
PTPRJ		PTPRJ	Confirmed by one method
Wee1		Wee1	
AKAP11			Confirmed by neither method
AKT2			
GSK3A			Confirmed by neither method
IRAK3			
ITGB1BP1			Confirmed by neither method
MAG3			
STK17A			Confirmed by neither method
STK39			
C6orf199		C6orf199	Confirmed by neither method
CSF1R			
ERBB4	ERBB4	ERBB4	Confirmed by neither method
ERN1			
FGFR2	FGFR2	FGFR2	Confirmed by neither method
K DR	K DR	K DR	
LMTK3	LMTK3	LMTK3	Confirmed by neither method
MAP2K7			
RAPGEF3			

High confidence hits → There is significant overlap between the hits that are confirmed by the two validation methods (≥ 2 siRNA and specificity-enhanced pools).


Low confidence hits (potential false positives) → Majority of hits with only one positive standard siRNA are not confirmed by specificity enhanced pools.

Low confidence hits (potential false positives) → A number of hits are either not expressed or are expressed at low levels.

Ambiguous hits (potential false positives/negatives) → Few hits are confirmed by only one method.

Seed controls for ruling out false negatives


A seed control is a chimeric siRNA that contains the 6 nucleotide seed sequence from the standard siRNA in a GAPDH siRNA scaffold that does not induce the phenotype (5).



IRAK3: Confirmed only by one standard siRNA; negative with the specificity-enhanced pool. The seed control is positive in phenotypic assay → the phenotype is due to seed-mediated off-targeting.


Specificity-enhancing chemical modification on the seed control abolishes the phenotype.

Low-confidence hits are due to seed-mediated off-target effects



Eight out of nine seed controls examined for low confidence hits induce seed-mediated phenotype. There was no seed-mediated effects on phenotype for the tested ambiguous hits. CENPE is a high confidence hit used as a control with no seed-mediated off-targeting.

Proposed hit stratification strategies



Follow up of candidate targets from primary screens with either deconvolution of standard pools or specificity-enhanced pools are equally viable approaches to prioritize potential targets based on confidence level.

High priority hits are identified by both methods
Low priority hits are negative by both methods
A few hits are ambiguous and confirmed by one or the other method

Seed-matched controls reveal that where only one siRNA produces the phenotype, often this is due to seed-mediated off targeting.

Seed-matched controls are important tools for ruling out seed-mediated off-target events during the hit stratification process.

References

1. A. Reynolds, et al., Rational siRNA design for RNA interference. *Nat. Biotechnol.*, 22(3), 326-30 (2004).
2. A.L. Jackson, et al., Position-specific chemical modification increases specificity of siRNA-mediated gene silencing. *RNA*, 12:71197-3205 (2006).
3. K.J. Simpson, et al., Identification of genes that regulate epithelial cell migration using an siRNA screening approach. *Nat. Cell. Biol.* 10, 1027-1038 (2008).
4. A. Birmingham, et al., 3' UTR seed matches, but not overall identity, are associated with RNAi off-targets. *Nat. Methods*, 3:3, 199-204 (2006).
5. E. Anderson, et al., Experimental validation of the importance of seed frequency to siRNA specificity. *RNA*, 14:5 (2008).

Acknowledgements

Special thanks to the R&D and production team at the Dharmacon site, now part of GE Healthcare, and the members of the RNAi Global Initiative.

gelifsciences.com/dharmacon

Dharmacon™
part of GE Healthcare

