


Dharmacon™ siRNA Selection Guide

We offer a wide selection of pre-designed siRNA product lines and formats. This quick selection guide will help to determine the best option for your particular needs.



Product Formats:

SMARTpool

A mixture of 4 siRNA provided as a single reagent; providing advantages in both potency and specificity.

Set of 4

A convenient option for purchasing aliquots of all 4 individual siRNAs targeting a single gene.

Set of 4 Upgrade

Discounted price on the Set of 4 based on a concurrent or prior purchase of the SMARTpool reagent for the same gene.

Individual siRNAs

Select 1, 2, 3 or 4 individual siRNAs per gene.

Dharmacon™ ON-TARGETplus™ and siGENOME™ SMARTpool™ reagents and 3 out of 4 siRNAs are guaranteed to silence mRNA by at least 75%. See website for details.

Ready to learn more about these offerings?

Contact your local Dharmacon representative or email us at ts.dharmacon@ge.com

gelifesciences.com/dharmacon