

YKO Parental Strains

Cat. #YSC1048, YSC1049, YSC1050, YSC1061, YSC1062

A nearly complete set of yeast open reading frame (ORF) knock-outs has been produced by the *Saccharomyces* Genome Deletion Project (SGDP) using five Yeast Knockout (YKO) Parental Strains derived from *Saccharomyces cerevisiae* S288C.¹ A PCR-based strategy was used to replace each ORF with a KanMX cassette containing unique tags, “barcodes”, for each deletion. Four different mutant collections have been generated: haploids of mating types MAT α and MAT α , homozygous diploids for non-essential genes, and heterozygous diploids, which contain the essential and non essential ORFs.

Storage:

- 4 °C for up to one week
- -80 °C indefinitely

Product Description:

Culture of *S. cerevisiae* in YPD broth with an inert growth indicator + 15% glycerol.

Genotypes of Specific Yeast Reference Strains¹:

BY4730	MAT α leu2 Δ 0 met15 Δ 0 ura3 Δ 0
BY4739	MAT α leu2 Δ 0 lys2 Δ 0 ura3 Δ 0
BY4741	MAT α his3 Δ 1 leu2 Δ 0 met15 Δ 0 ura3 Δ 0
BY4742	MAT α his3 Δ 1 leu2 Δ 0 lys2 Δ 0 ura3 Δ 0
BY4743	4741/4742

Homozygous diploids are in the BY4743 background unless 4730/4739 is indicated.

Useful Websites and References:

Saccharomyces Genome Deletion Project, click [here](#).

1. Winzeler, E.A., et al., Functional characterization of the *Saccharomyces cerevisiae* Genome by gene deletion and parallel analysis, *Science*, 285:901-906, 1999.

GE, imagination at work and GE monogram are trademarks of General Electric Company. Dharmacon is a trademark GE Healthcare companies. All other trademarks are the property of General Electric Company or one of its subsidiaries. ©2014 General Electric Company—All rights reserved. Version published June 2014. GE Healthcare UK Limited, Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

